

Iron and Steel

Ferrous Salt in Etching Solution

Redox titration by
Automatic Potentiometric Titrator

Standard

1. Abstract

Here we demonstrate measurement of concentration of ferrous salt in etching solution according to the method from technical literatures as referred below.

When a sample does not contain chlorine ion, dissolve it in sulfuric solution, and titrate with 0.01mol/L potassium permanganate.



If it contains chlorine ion, add manganese sulfate to void side reaction. If it contains too much hydrochloric acid, use vacuum evaporative condenser to dispel it.

2. Reference

- 1) Experiment and Calculation in Quantitative Analysis - Vol.2 - Kyoritsu Publication

3. Cautions in measurement

- 1) When ferrous salt in hydrochloric acid solution is titrated with potassium permanganate, the potassium permanganate is consumed as per below formula, generating chlorine or hypochlorous acid, which causes measurement error. It is assumed the generation of chlorine is due to hydrochloric acid strongly oxidized by manganese ions which exist temporarily during chemical reaction.



Chlorine is generated even in dilute hydrochloric acid if ferric ion coexists. This is supposedly by concurrent reaction of ferric ions. However, if manganese ion coexists, chlorine is not generated. It is assumed that highly oxidizing manganese ions oxidize ions II to III and IV, then ferric ions, and returns to original manganese ions before oxidizing hydrochloric acid and that manganese ion II reduces redox potential of potassium permanganates. It is practiced in the analysis of irons ore to titrate with potassium permanganate by intentionally adding manganese ions II while maintaining hydrochloric acidity. In this test, the sample in etching solution contains chlorine ions, and manganese sulfate is added in order to avoid side reaction.

- 2) For manual analysis, titrate with potassium permanganate until it shows reddish color. Dense sample may show yellowish color due to ferrous ions generated in titration, however, it turns to colorless when phosphoric acid is added, thus generating complex ion of $\text{Fe}(\text{PO}_4)_3^{6-}$.

4. Post-measurement care

Rinse the electrode with water, and keep it dipped in pure water to avoid drying up.

5. Test equipment

Main unit : Automatic potentiometric titrator (Standard preamplifier: STD—)

Electrode : Option Combination Pt. electrode

6. Reagent

Reagent : 0.01mol/L(0.05N) Potassium permanganate solution

Additive : Manganese sulfate (Crystalline manganese sulfate, Phosphoric acid, Sulfuric acid)

7. Measurement procedure

—Preparation—

- 1) Add 200mL pure water to 90g crystalline manganese sulfate ($\text{MnSO}_4 \cdot 4\text{H}_2\text{O}$) to dissolve, and add 175mL phosphorous acid (SG:1.7), then add 175mL SO_4 (SG:1.84) gradually, and cool it. Add water to make it total 1L, and leave it for 3 days. The above is how the manganese sulfate was prepared.

—Measurement—

- 1) Weigh 0.5g sample and put it in a 200mL beaker.
- 2) Add 40mL manganese sulfate and 50mL pure water.
- 3) Titrate with 0.01mol/L potassium permanganate to obtain concentration of ferrous salt (FeCl_2).
- 4) Perform a blank test under the same conditions to correct titration volume.

8. Formula

Concentration of ferrous salt (%) = $(\text{EPl} - \text{Bl1}) \times \text{TF} \times \text{C1} \times \text{K1} / \text{SIZE}$

EPl : Titration volume (mL)

Bl1 : Blank level (0.0108mL)

TF : Factor of reagent (1.008)

C1 : Concentration conversion coefficient (6.3377mg/mL)
(FeCl_2 equivalent to 1mL of 0.01mol/L KMnO_4)

K1 : Unit conversion coefficient (0.1)

SIZE : Sample size (g)

9. Example of measurement

— Ambient condition —

Room temperature : 23 °C	Humidity : 45 %	Weather : Fair
--------------------------	-----------------	----------------

-Titration parameters-

Model : AT-510	
Method No. : 17	
Titr.mode : Auto	
Intermit	
Titr.form : EP Stop	
[Titration parameter]	[Result parameter]
Form : EP Stop	<Calculation>
APB No. : 1	Calc.Type : Sample
Unit No. : 1	Conc.1 : Set
Detector No. : 1	CO1=
Unit : mV	(EP1-BL1)*TF*C1*K1/SIZ
Max. Volume : 20.0mL	E
Wait Time : 0s	Unit : %
Direction : Auto	EP No. : 1
	Temp.Comp. : Off
[Control parameter]	<Constant>
End Point No. : 1	C1 : 6.3377
End Sense : Auto	K1 : 0.1
End Point Area : Off	<Titr. Constant>
Separation : Off	Factor : 1.008
Over Titr.Vol. : 0mL	<Blank>
Gain : 1	Blank1 : 0.0108
Data samp.Pot. : 4.0mV	
Data samp.Vol. : 0.5mL	
Stability : 0.5mV/s	
Delay Time : 1s	
Limit Time : 30s	

-Titration curve-

*** Result ***
Sample No. : 03-01
Date : 2001/06/04 16:27
Sample ID :
Method No. : 17
<Auto Intermit>
Method Name :
Titr.time : 00:11:19
Size : <u>0.5163g</u>
Conc-1 : <u>13.108%</u>
End point-1
Volume : <u>10.6041mL</u>
Potential : 680.0mV

(The above printout data were obtained from titration by AT-510 unit.)

«Titration parameter»

Form: of formula / APB No. the burette used in titration / Unit No. [APB Unit File number](#)
 Detector No. the detector used in titration / Unit: potential unit in EP detection / Max Volume: of titration
 Wait Time: before titration starts / Direction: of titration

«Control parameter»

End Point No.: number of EPs / End sense: of endpoint / End Point Area: detection area / Separation: of potential
 Over Titr.Vol.: over-titration / Gain: sensitivity of detection signal / Data samp.Pot.: potential changes of sampling signal
 Data samp.Vol. titration volume for sampling detection signal / Control Speed: of dosing / Stability: of EP sense
 Delay Time: before stability check / Limit Time: for stability check

«Result parameter»

Calc.Type of formula / Conc.1: formula 1 / Unit: of result / EP No. of calculation / Temp.Comp.: temperature compensation
 C1: concentration coefficient / K1: unit conversion coefficient / Factor: of reagent / Blank1: blank level 1

–Measurement results–

N	Sample (g)	Titration (mL)	FeCl ₂ (%)
1	<u>0.5163</u>	<u>10.6041</u>	<u>13.108</u>
2	0.5199	10.4438	12.820
3	0.5096	10.5516	13.214

Concentration of FeCl ₂	
Mean	13.047 %
SD	0.2039 %
RSD	1.5627 %

* The above results were obtained by 3 tests of the same sample.

* Red underline shows the data from page 3/4.

10. Summary

The sample measurement shows favorable results of a good repeatability with 1.5% relative standard deviation.

Precise and reliable measurement is assured by the automated potentiometry.