

Petroleum

## Density of Low Sulphur Waxy Residue (tar oil)

Density measurement by  
Density/Specific Gravity Meter

Standard	JIS K 0061	ASTM D 4052	API Std. 2540
	JIS K 2249	ASTM D 5002	
	ASTM D1250	ISO 91-1	

### 1. Abstract

Measurement of density of crude petroleum and petroleum products is specified in JIS K 2249-1995 Crude petroleum and petroleum products - Determination of density and petroleum measurement tables based on a reference temperature.

Density of crude petroleum and petroleum products is prescribed to express value [density (15°C)] in 15°C normally. The sample which is hard to measure such as solid at room temperature or high viscosity liquid is measured after lowering viscosity by heating with a method converted to density at room temperature. This application note exemplifies a measurement of specific gravity 60/60°F and specific gravity 15/15°C with a built-in temperature conversion formula using multiple sample changer (high temperature specification (electric heater type)) and density/specific gravity meter made by Kyoto Electronics (KEM).

This density/specific gravity meter is equipped conversion formula to density, API degree and specific gravity (temperature: 60°F•15°C•20°C) of crude petroleum(A)•fuel oil(B)•grease(D) based on prescription of petroleum products.

### 2. Reference

- 1) JIS K 0061-2001 Test methods for density and relative density of chemical products
- 2) JIS K 2249-1995 Crude petroleum and petroleum products - Determination of density and petroleum measurement tables based on a reference temperature
- 3) ASTM D1250-08 Standard Guide for Use of the Petroleum Measurement Tables
- 4) ASTM D 4052-09 Standard Test Method for Density, Relative Density, and API Gravity of Liquids by Digital Density Meter
- 5) ASTM D 5002-99(2010) Standard Test Method for Density and Relative Density of Crude Oils by Digital Density Analyzer
- 6) ISO 91-1:1992 Petroleum measurement tables -- Part 1: Tables based on reference temperatures of 15 degrees C and 60 degrees F
- 7) API Std. 2540

### 3. Cautions in measurement

- 1) Use a fresh desiccant (silica gel). If its color is reddish, change it with new one (blue).
- 2) Before sample measurement, perform factor calibration with dry air and pure water after degas.
- 3) Watch out for fire due to use flammable solvent.

## 4. Post-measurement care

Drain out residue inside the cell. Rinse well the measurement cell with solvent, and dry it completely.

## 5. Test equipment

Main unit : Density/specific gravity meter

Option : Multiple sample changer (high temperature specification (electric heater type))

## 6. Reagent

Rinse solvent : Toluene (for cleaning) rinse solvent 1

Rinse solvent : Acetone (for drying) rinse solvent 2

## 7. Measurement procedure

—Calibration—

- 1) Set the temperature of density meter to 70°C, and leave it until stabilized.
- 2) Set the temperature of turn table side of multiple sample changer to 80°C which is 10°C higher than measurement temperature and tube heating unit to 75°C which is 5°C higher than measurement temperature, and leave it until stabilized.
- 3) Calibrate the meter with dry air and pure water after degas. The water should be used after heated at 100 °C.

—Measurement—

- 1) Set sample to sampling bottle and multiple sample changer, and make sampling temperature stable.
- 2) Start measurement.  
(when measurement is started, the series of sampling, measurement, cleaning and drying are automatically operated.)

## 8. Example of measurement

– Ambient condition –

Room temperature : 25 °C	Humidity : 54 %	Weather : Fair
--------------------------	-----------------	----------------

<b>-Measurement Parameter-</b>	
[Method Parameter]	
<Measurement Parameter>	
Temperature	: 70.00 °C
Stability	: 1
Limit Time	: 600 s
Visco. Corr.	: Off
Sequence	: Off
Calib.Mate.	: Air&Water
Sample[mL]	: 0.00
<Contents>	
Contents Name	: Density
Decimal	: 4
Contents Name	: API(S.G.)A 60F API(S.G.)A 15C
Decimal	: 4
<Temperature Comp.>	
Temp.Comp.	: Off
<Sequence>	
Sequence File	: Sequence auto No.
01	Sampling
02	Meas.
03	Drain
04	Rince1
05	Rince2
06	Purge
Sampling Seq.	: Auto
Samp.Limit	: 0s
O.S.Rate	: 70%
Drain Seq.	: Auto
Drain Rate	: 100%
Rince-1 Time	: 30s
Rince-2 Time	: 10s
Purge Seq.	: Auto
Tolerance	: 10
<b>-Calibration Parameter-</b>	
[Check&Calib]	
One Point Calib.	: off
Calib. Temp.	: 70.00°C
Stability	: 1
Viscosity	: Off
Limit Time	: 600s

(Rest of Calibration parameter)	
Tolerance	: 0.0002
Sequence	: On
Sequence Name	: Sequence 1
Calib. Mate.	: Air&Water
Check	: Off
<b>-Sequence-</b>	
<Sequence>	
Sequence File	: Sequence 1 No.
01	Sampling
02	Meas.
03	Drain
04	Rince2
05	Purge
Sampling Seq.	: Auto
O.S.Rate	: 70%
Samp.Limit	: 20s
Drain Seq.	: Auto
Drain Rate	: 100%
Rince-1 Time	: 0s
Rince-2 Time	: 10s
Purge Seq.	: Auto
Tolerance	: 10
<b>-Measurement Result-</b>	
*** R e s u l t ***	
Sample No.	01-001
Date	: 2011/04/12 15:41
Sample ID	:
Method Name	: original
Meas.Temp.	: 70.00 °C
d[g/cm3]	: <u>1.0577</u>
<Result>	
API(S.G.)A 60F	: <u>1.0897</u>
API(S.G.)A 15C	: <u>1.0900</u>
Meas.Time	: 00:01 : 19

(Printout examples by DA-640)

—Measurement result—

	Frequency	Density (g/cm <sup>3</sup> )	API(S.G.)A 60F	API(S.G.)A 15C	Statistics	
					Mean	1.0577 g/cm <sup>3</sup>
1	<u>112121</u>	<u>1.0577</u>	<u>1.0897</u>	<u>1.0900</u>	Diff (Max-Min)	0.0000 g/cm <sup>3</sup>
2	112121	1.0577	1.0897	1.0900	SD	0.0000 g/cm <sup>3</sup>
3	112122	1.0577	1.0897	1.0900	RSD	0.00 %
4	112123	1.0577	1.0897	1.0900		
5	112123	1.0577	1.0897	1.0900		

- \* The above data are the results of 5 tests of the same sample.
- \* Red underline shows the data from page 3/4.
- \* The above "API(S.G.)A 60F" shows gravity 60/60°F.
- \* The above "API(S.G.)A 15C" shows gravity 15/15°C.

## 9. Summary

Sample measurement this time shows a good repeatability.

Since the following density/specific gravity meters are equipped with density conversion formula corresponding to JIS K2249·ASTM D1250·ASTM D 4052·ASTM D 5002·ISO91·API Std.2540, can be calculated gravity 60/60°F and gravity 15/15°C of crude petroleum(A)·fuel oil(B)·grease(D) automatically.